



ULTRASOUND    PROBE REPAIR



We are an internationally active company with a focus on ultrasonography. Priority is given to the provision of tailored solutions for clinical operations and specialist physician practices.

Expert advice, high-quality products and services and needs-based support create a high degree of customer satisfaction and thereby contribute to shared success.

An in-house research and development department devoted to probe repairs make us a long-term and quality-oriented partner in the medical technology sector.



## PRODUCTS AND SERVICES

### **Ultrasound systems:**

Advisory service, sales and distribution and servicing of ultrasound systems  
New, demo or used


### **Probe repair:**

Professional repair of all types of ultrasound probes from all manufacturers

### **Refurbishment:**

Quality assurance measures for your ultrasound systems/probes.

# PROBE ANALYSIS



A comprehensive initial analysis forms the basis for a successful repair of your probes. Timely investigation can discover defects before these necessitate expensive repairs. Our service offering in the area of probe analysis comprises:

- Check-in, optical and technical
- Advance detailed probe analysis
- Quality-assured probe repair
- Post-repair detailed probe analysis
- Test protocol preparation
- Check-out, optical and technical
- Service warranty
- Loan probes

## PROBE REPAIR



**Standard probes** (curved, linear and phased array probes):

- Lens coating renewal
- Cable anti-kink protector renewal
- Probe housing processing and repair
- Connector repair
- Cable replacement and repair



**TEE probes:**

- IT-coating renewal/replacement
- Bending rubber replacement
- Pulley renewal and calibration
- Lens replacement



**3D volume probes:**

- Lens cap replacement and renewal
- Probe mechanics cleaning and processing
- Probe housing processing and repair
- Cable replacement and repair



RESEARCH • DEVELOPMENT • INNOVATION



Our repair laboratory continuously develops new procedures for the optimisation of probe repair.

FAST • COMPETENT • DIRECT

Nationwide on-site support provided by our trained specialist staff.



EXPECT • MORE

Ortner Medical Solutions

Gewerbeparkstraße 3 • A-8143 Dobl bei Graz

M: office@ort-med.com • T: +43 3136 62155 • F: +43 3136 62155 4

[www.ort-med.com](http://www.ort-med.com)